

ALDH7A1 antibody



Catalog Number: 107974

Product name

ALDH7A1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ALDH7A1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in mouse liver tissue, HepG2 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of ALDH7A1 recombinant protein (Accession Number: NM_001201377). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107974(ALDH7A1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:107974(ALDH7A1 antibody) at dilution of 1:50 (under 40x lens)