

# ALK,CD246 antibody



Catalog Number: 107983

## Product name

ALK,CD246 antibody

## Immunogen

[Human ALK,CD246 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ALK,CD246 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human brain tissue.

## Preparation

This antibody was obtained by immunization of ALK,CD246 recombinant protein (Accession Number: NM\_004304). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IHC


## Dilutions

Recommended Dilution:  
IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:107983(ALK,CD246 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:107983(ALK,CD246 Antibody) at dilution of 1:50 (under 40x lens)