

AMIGO2 antibody



Catalog Number: 107993

Product name

AMIGO2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AMIGO2 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, K-562 cells. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of AMIGO2 recombinant protein (Accession Number: XM_047428785). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:107993(AMIGO2 antibody) at dilution of 1:500