

# AMPK beta 2 antibody



Catalog Number: 108009

## Product name

AMPK beta 2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

AMPK beta 2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, fetal human brain tissue, HeLa cells, human heart tissue, mouse testis tissue. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of AMPK beta 2 recombinant protein (Accession Number: NM\_005399). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108009(PRKAB2 antibody) at dilution of 1:400