

**Product name**

Amyloid beta (A4) precursor protein antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

Amyloid beta (A4) precursor protein Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 102 kDa

**Preparation**

This antibody was obtained by immunization of Amyloid beta (A4) precursor protein recombinant protein (Accession Number: BC004369). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG

**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Validations**

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108016(APP antibody) at dilution of 1:400