

# ANKZF1 antibody



Catalog Number: 108087

## Product name

ANKZF1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ANKZF1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, COLO 320 cells, HeLa cells, HepG2 cells, MCF-7 cells. Observed molecular weight by Western-blot: 81 kDa

## Preparation

This antibody was obtained by immunization of ANKZF1 recombinant protein (Accession Number: NM\_001042410). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108087(ANKZF1 Antibody) at dilution of 1:2000