

ANO8 antibody

Catalog Number: 108103

Product name

ANO8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ANO8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 136kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020959). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108103(ANO8 antibody) at dilution of 1:500



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:108103(ANO8 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG