

# AP1,JUN,P39 antibody



Catalog Number: 108114

## Product name

AP1,JUN,P39 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AP1,JUN,P39 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells, rat brain tissue. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 40-46 kDa

## Preparation

This antibody was obtained by immunization of AP1,JUN,P39 recombinant protein (Accession Number: BC002646). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108114(AP1,JUN,P39 Antibody) at dilution of



Immunofluorescent analysis of (10% Formaldehyde) fixed NIH/3T3 cells using Catalog No:108114(AP1,JUN,P39 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)