

AP2A2 antibody

Catalog Number: 108119

Product name

AP2A2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP2A2 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 104 kDa

Preparation

This antibody was obtained by immunization of AP2A2 recombinant protein (Accession Number: NM_012305). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108119(AP2A2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108119(AP2A2 Antibody) at dilution of 1:50 (under 40x lens)