

ARMC5 antibody



Catalog Number: 108150

Product name

ARMC5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ARMC5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A2780 cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of ARMC5 recombinant protein (Accession Number: BC014945). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108150(ARMC5 Antibody) at dilution of 1:300



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108150(ARMC5 Antibody) at dilution of 1:300