

Product name

ARPC1B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ARPC1B Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer.

Preparation

This antibody was obtained by immunization of ARPC1B recombinant protein (Accession Number: XM_006715825). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:108164(ARPC1B antibody) at dilution of 1:100 (under 25x lens)