

ARHGAP10 antibody



Catalog Number: 108260

Product name

ARHGAP10 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ARHGAP10 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Positive WB detected in human heart tissue, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 89 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_024605). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108260(ARHGAP10 antibody) at dilution of



1:800
Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:108260(ARHGAP10 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:108260(ARHGAP10 antibody) at dilution of

ARHGAP10 antibody



Catalog Number: 108260

1:50 (under 40x lens)