

# ATG13 antibody



Catalog Number: 108290

## Product name

ATG13 antibody

WB: 1:500-1:5000

IP: 1:200-1:1000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

ATG13 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in human brain tissue, BGC-823 cells, mouse cerebellum tissue, mouse thymus tissue, SH-SY5Y cells. Positive IP detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 57-63kd

## Validations




## Preparation

This antibody was obtained by immunization of ATG13 recombinant protein (Accession Number: NM\_001346342). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108290(KIAA0652 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human brain using Catalog No:108290(KIAA0652 antibody) at dilution of 1:100 (under 10x lens)

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:



# ATG13 antibody

Catalog Number: 108290

---

Immunohistochemical of paraffin-embedded human brain using Catalog No:108290(KIAA0652 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-KIAA0652 (IP:Catalog No:108290, 4ug; Detection:Catalog No:108290 1:300) with SH-SY5Y cells lysate 1600ug.

