

ATRIP antibody

Catalog Number: 108327

Product name

ATRIP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATRIP Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, A549 cells, BxPC-3 cells, HeLa cells, NIH/3T3 cells, PC-3 cells, RAW264.7 cells. Positive IP detected in MCF-7 cells.

Observed molecular weight by Western-blot: 86-90kd

Preparation

This antibody was obtained by immunization of ATRIP recombinant protein (Accession Number: NM_032166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108327(ATRIP antibody) at dilution of 1:1000



IP Result of anti-ATRIP (IP:Catalog No:108327, 3ug; Detection:Catalog No:108327 1:1000) with MCF-7 cells lysate 2800ug.