

ATP4B antibody



Catalog Number: 108349

Product name

ATP4B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ATP4B Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, rat kidney tissue. Observed molecular weight by Western-blot: 60-80 kDa

Preparation

This antibody was obtained by immunization of ATP4B recombinant protein (Accession Number: NM_000705). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:108349(ATP4B antibody) at dilution of 1:800