

# B4GALT7 antibody



Catalog Number: 108397

---

## Product name

B4GALT7 antibody

## Specificity

human; other species not tested.

## Antibody description

B4GALT7 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells. Observed molecular weight by Western-blot: 37kd

## Preparation

This antibody was obtained by immunization of B4GALT7 recombinant protein (Accession Number: NM\_007255). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108397(B4GALT7 antibody) at dilution of 1:500