

# BEND4 antibody

Catalog Number: 108444

## Product name

BEND4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BEND4 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human heart tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 55-58 kDa; 33-36 kDa

## Preparation

This antibody was obtained by immunization of BEND4 recombinant protein (Accession Number: BC128592). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117117(BEND4 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:117117(BEND4 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:117117(BEND4 Antibody) at dilution of 1:50

# BEND4 antibody



Catalog Number: 108444

---

(under 40x lens)