

BLOC1S3 antibody



Catalog Number: 108454

Product name

BLOC1S3 antibody

Specificity

Human; other species not tested.

Antibody description

BLOC1S3 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human cervical cancer tissue.

Preparation

This antibody was obtained by immunization of BL°C1S3 recombinant protein (Accession Number: NM_212550). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:117189(BLoc1S3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:117189(BLoc1S3 Antibody) at dilution of 1:50 (under 40x lens)