

# BMPR2 antibody

Catalog Number: 108468

## Product name

BMPR2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BMPR2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 120kd, 97kd

## Preparation

This antibody was obtained by immunization of BMPR2 recombinant protein (Accession Number: NM\_001204). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117203(BMPR2 antibody) at dilution of 1:1000



IP Result of anti-BMPR2 (IP:Catalog No:117203, 4ug; Detection:Catalog No:117203 1:500) with mouse brain tissue lysate 6500ug.