

BMPR2 antibody



Catalog Number: 108468

Product name

BMPR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BMPR2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 120kd, 97kd

Preparation

This antibody was obtained by immunization of BMPR2 recombinant protein (Accession Number: NM_001204). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117203(BMPR2 antibody) at dilution of 1:1000



IP Result of anti-BMPR2 (IP:Catalog No:117203, 4ug; Detection:Catalog No:117203 1:500) with mouse brain tissue lysate 6500ug.