

BOLA2 antibody



Catalog Number: 108477

Product name

Polyclonal

BOLA2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

BOLA2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human testis tissue. Observed molecular weight by Western-blot: 7 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of BOLA2 recombinant protein (Accession Number: NM_001031833). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117212(BOLA2 Antibody) at dilution of 1:1000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality