

# AXIN1 antibody



Catalog Number: 108517

---

## Product name

AXIN1 antibody

## Immunogen

[Human AXIN1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

AXIN1 Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells, A431 cells, A549 cells. Observed molecular weight by Western-blot: 100-110kd

## Preparation

This antibody was obtained by immunization of AXIN1 recombinant protein (Accession Number: NM\_181050). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:117078(AXIN1 antibody) at dilution of 1:500