

AXIN2 antibody



Catalog Number: 108518

General Information

Product name:

AXIN2 antibody

Immunogen:

[Human AXIN2 Recombinant protein \(GST tag\)](#)

Applications:

ELISA, WB

Specificity:

human, mouse, rat; other species not tested.

Antibody description:

AXIN2 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse lung tissue. Observed molecular weight by Western-blot: 93-100kd

Properties

Preparation:

This antibody was obtained by immunization of AXIN2 recombinant protein (Accession Number: BC006295). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

Recommended Dilution:

WB: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:117079(AXIN2 antibody) at dilution of 1:200