

# B3GALT4 antibody



Catalog Number: 108527

## Product name

B3GALT4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

B3GALT4 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human pancreas cancer tissue.

## Preparation

This antibody was obtained by immunization of B3GALT4 recombinant protein (Accession Number: NM\_003782). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:117088(B3GALT4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:117088(B3GALT4 antibody) at dilution of 1:100 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117088(B3GALT4 Antibody) at dilution of 1:1000