

B3GNT3 antibody



Catalog Number: 108530

Product name

B3GNT3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

B3GNT3 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, mouse lung tissue. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 43kd

Preparation

This antibody was obtained by immunization of B3GNT3 recombinant protein (Accession Number: XM_011527626). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:117091(B3GNT3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon using Catalog No:117091(B3GNT3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:117091(B3GNT3 antibody) at dilution of 1:100 (under 40x lens)