

BTBD8 antibody



Catalog Number: 108538

Product name

BTBD8 antibody

Specificity

Human; other species not tested.

Antibody description

BTBD8 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 43, 34 kDa

Preparation

This antibody was obtained by immunization of BTBD8 recombinant protein (Accession Number: NM_183242). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108538(BTBD8 Antibody) at dilution of 1:300