

# BTNL2 antibody

Catalog Number: 108549

---

## Product name

BTNL2 antibody

## Immunogen

[Human BTNL2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BTNL2 Rabbit Polyclonal antibody. Positive WB detected in mouse small intestine tissue, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of BTNL2 recombinant protein (Accession Number: BC127642). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:108549(BTNL2 Antibody) at dilution of 1:600