

BUB3 antibody

Catalog Number: 108552

Product name

BUB3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

BUB3 Rabbit Polyclonal antibody. Positive IF detected in HepG2.

Preparation

This antibody was obtained by immunization of BUB3 recombinant protein (Accession Number: NM_004725). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of HepG2 cells, using BUB3 antibody Catalog No:108552 at 1:100 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent DNA dye).