

# C10orf27 antibody



Catalog Number: 108562

## Product name

C10orf27 antibody

## Immunogen

[Human C10orf27 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

C10orf27 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 39 kDa

## Preparation

This antibody was obtained by immunization of C10orf27 recombinant protein (Accession Number: NM\_001318242). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108562(C10orf27 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:108562(C10orf27 Antibody) at dilution of 1:50