

BRF2 antibody



Catalog Number: 108587

Product name

BRF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BRF2 Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells, DU 145 cells. Observed molecular weight by Western-blot: 50-60 kDa

Preparation

This antibody was obtained by immunization of BRF2 recombinant protein (Accession Number: NM_018310). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:117234(BRF2 antibody) at dilution of 1:300