

# BRMS1L antibody



Catalog Number: 108593

## Product name

BRMS1L antibody

## Specificity

Human; other species not tested.

## Antibody description

BRMS1L Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, Hela cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of BRMS1L recombinant protein (Accession Number: NM\_032352). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117240(BRMS1L Antibody) at dilution of 1:600