

C16orf14 antibody



Catalog Number: 108604

Product name

C16orf14 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C16orf14 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Positive IF detected in Hela cells, HepG2 cells. Observed molecular weight by Western-blot: 8-10 kDa

Preparation

This antibody was obtained by immunization of C16orf14 recombinant protein (Accession Number: BC007346). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:108604(C16orf14 antibody) at dilution of



1:2000

Immunofluorescent analysis of Hela cells, using C16orf14 antibody Catalog No:108604 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).