

C16orf70 antibody



Catalog Number: 108608

Product name

C16orf70 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C16orf70 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, L02 cells. Positive IHC detected in human liver cancer tissue, human lung tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of C16orf70 recombinant protein (Accession Number: NM_001320541). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108608(C16orf70 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:108608(C16orf70 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:108608(C16orf70 antibody) at dilution of 1:100

C16orf70 antibody



Catalog Number: 108608

(under 40x lens)