

C16orf93 antibody

Catalog Number: 108611

Product name

C16orf93 antibody

Immunogen

[Human C16orf93 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C16orf93 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of C16orf93 recombinant protein (Accession Number: NM_001014979). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108611(C16orf93 Antibody) at dilution of 1:700