

# C18orf45 antibody



Catalog Number: 108615

## Product name

C18orf45 antibody

## Specificity

Human; other species not tested.

## Antibody description

C18orf45 Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 35-40 kDa

## Preparation

This antibody was obtained by immunization of C18orf45 recombinant protein (Accession Number: NM\_032933). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:108615(C18orf45 Antibody) at dilution of 1:600