

C19orf34 antibody



Catalog Number: 108620

Product name

C19orf34 antibody

Specificity

Human; other species not tested.

Antibody description

C19orf34 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of C19orf34 recombinant protein (Accession Number: NM_152771). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108620(C19orf34 Antibody) at dilution of 1:500