

# C10orf67 antibody



Catalog Number: 108634

## Product name

C10orf67 antibody

## Specificity

Human; other species not tested.

## Antibody description

C10orf67 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, A549 cells, human testis tissue. Observed molecular weight by Western-blot: 18-21 kDa

## Preparation

This antibody was obtained by immunization of C10orf67 recombinant protein (Accession Number: NM\_153714). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108634(C10orf67 Antibody) at dilution of 1:300