

C15orf23 antibody



Catalog Number: 108666

Product name

C15orf23 antibody

Specificity

Human; other species not tested.

Antibody description

C15orf23 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HL-60 cells. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of C15orf23 recombinant protein (Accession Number: NM_033286). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108666(C15orf23 Antibody) at dilution of 1:1500