

# C15orf39 antibody



Catalog Number: 108669

## Product name

C15orf39 antibody

## Specificity

Human; other species not tested.

## Antibody description

C15orf39 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A549 cells, HeLa cells. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of C15orf39 recombinant protein (Accession Number: BC001762). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108669(C15orf39 Antibody) at dilution of 1:600