

C1orf226 antibody

Catalog Number: 108674

Product name

C1orf226 antibody

Specificity

Human; other species not tested.

Antibody description

C1orf226 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of C1orf226 recombinant protein (Accession Number: NM_001085375). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-C1orf226 (IP:Catalog No:108674, 3ug; Detection:Catalog No:108674 1:300) with A431 cells lysate 1520ug.