

C1orf77 antibody



Catalog Number: 108681

Product name

C1orf77 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

C1orf77 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 26-30 kDa

Preparation

This antibody was obtained by immunization of C1orf77 recombinant protein (Accession Number: BC002733). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108681(C1orf77 antibody) at dilution of 1:500