

C3 / C3b / C3c antibody



Catalog Number: 108718

Product name

C3 / C3b / C3c antibody

IP: 1:500-1:5000

IHC: 1:50-1:500

Specificity

Human; other species not tested.

Antibody description

C3 / C3b / C3c Rabbit Polyclonal antibody. Positive IP detected in human plasma tissue. Positive WB detected in human blood tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 115 kDa

Preparation

This antibody was obtained by immunization of C3 / C3b / C3c recombinant protein (Accession Number: NM_000064). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:108718(C3 antibody) at dilution of 1:1000



IP Result of anti-C3 (IP:Catalog No:108718, 4ug; Detection:Catalog No:108718 1:1000) with human plasma tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog

C3 / C3b / C3c antibody

Catalog Number: 108718

No:108718(C3 Antibody) at dilution of 1:100
(under 10x lens).



Immunohistochemistry of paraffin-embedded
human liver tissue slide using Catalog
No:108718(C3 Antibody) at dilution of 1:100
(under 40x lens).