

# C3orf20 antibody

Catalog Number: 108720

## Product name

C3orf20 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

C3orf20 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 101 kDa, 88 kDa

## Preparation

This antibody was obtained by immunization of C3orf20 recombinant protein (Accession Number: NM\_032137). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108720(C3orf20 antibody) at dilution of 1:500



IP Result of anti-C3orf20 (IP:Catalog No:108720, 3ug; Detection:Catalog No:108720 1:500) with mouse brain tissue lysate 10000ug.