

C3orf27 antibody

Catalog Number: 108722

Product name

C3orf27 antibody

Specificity

Human; other species not tested.

Antibody description

C3orf27 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of C3orf27 recombinant protein (Accession Number: NM_007354). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-C3orf27 (IP:Catalog No:108722, 4ug; Detection:Catalog No:108722 1:500) with HeLa cells lysate 1200ug.