

C4BPB antibody



Catalog Number: 108729

Product name

C4BPB antibody

Immunogen

Specificity

Human, Mouse; other species not tested.

Antibody description

C4BPB Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue, HeLa cells, mouse ovary tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of C4BPB recombinant protein (Accession Number: NM_001017364). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:108729(C4BPB Antibody) at dilution of 1:600