

C5orf24 antibody

Catalog Number: 108735

Product name

C5orf24 antibody

Specificity

Human; other species not tested.

Antibody description

C5orf24 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human kidney tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of C5orf24 recombinant protein (Accession Number: XM_017009050). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108735(C5orf24 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog No:108735(C5orf24 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog

C5orf24 antibody



Catalog Number: 108735

No:108735(C5orf24 Antibody) at dilution of 1:50 (under 40x lens)