

# Cadherin-18 antibody



Catalog Number: 108775

## Product name

Cadherin-18 antibody

## Specificity

Human; other species not tested.

## Antibody description

Cadherin-18 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, A549 cells. Observed molecular weight by Western-blot: 105kd

## Preparation

This antibody was obtained by immunization of Cadherin-18 recombinant protein (Accession Number: XM\_017008925). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:108775(CDH18 antibody) at dilution of 1:300