

C6orf204 antibody

Catalog Number: 108783

Product name

C6orf204 antibody

Specificity

Human; other species not tested.

Antibody description

C6orf204 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human tonsillitis tissue. Positive WB detected in PC-13 cells. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of C6orf204 recombinant protein (Accession Number: NM_206921). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:108783(C6orf204 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108783(C6orf204 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108783(C6orf204 Antibody) at dilution of 1:50

C6orf204 antibody



Catalog Number: 108783

(under 40x lens)