

# C9orf41 antibody



Catalog Number: 108799

## Product name

C9orf41 antibody

## Immunogen

[Human C9orf41 Recombinant protein \(His tag\)](#)

## Specificity

human; other species not tested.

## Antibody description

C9orf41 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A549 cells, HeLa cells, MCF7 cells.

## Preparation

This antibody was obtained by immunization of C9orf41 recombinant protein (Accession Number: NM\_152420). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108799(C9orf41 antibody) at dilution of 1:1000