

Cadm1 antibody



Catalog Number: 108812

Product name

Cadm1 antibody

Specificity

mouse; other species not tested.

Antibody description

Cadm1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of Cadm1 recombinant protein (Accession Number: AF539424). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108812(Cadm1 antibody) at dilution of 1:500